

TECHNICAL BULLETIN

NUMBER 2021 - 004

Date: March 26, 2021
Author: Justin Dickerman
Subject: PeCOD Inorganic Interference Review

The following tables summarize the impact of a range of common inorganic anions and cations on the determination of COD using the PeCOD® technique.

For each inorganic species, solutions containing 0, 20, 50, 100, 250, 500 ppm (by mass) of the anion or cation, 60ppm COD (as sorbitol) and 1M LiNO₃ (containing 20ppm COD spike) were prepared and analyzed, unless otherwise stated. Therefore, the below ion concentrations represent the concentration IN THE CELL (i.e. if analyzed in a different range, the interference levels may vary due to different electrolyte dilution effects).

Anions	Formula	Remark
Ammonium	NH ₄ ⁺	No interference for NH ₄ ⁺ ≤ 500 ppm Note: Similar results are obtained for Ammonia
Carbonate	CO ₃ ²⁻	No interference for CO ₃ ²⁻ ≤ 500ppm using chloride resistant sensor
Chlorate	ClO ₃ ⁻	No interference for ClO ₃ ⁻ ≤ 500ppm
Chloride	Cl ⁻	No interference for Cl ⁻ < 200 ppm. COD reduced by up to 20% at Cl ⁻ levels of 500ppm using Chloride resistant sensor. Other halides (F-, Br-, I-) would be expected to behave in the same manner
Nitrate	NO ₃ ⁻	No interference, NO ₃ ⁻ can be used as PeCOD electrolyte
Nitrite	NO ₂ ⁻	No interference for NO ₂ ⁻ ≤ 500ppm
Perchlorate	ClO ₄ ⁻	No interference, ClO ₄ ⁻ can be used as PeCOD electrolyte
Phosphate	PO ₄ ³⁻	No interference for PO ₄ ³⁻ ≤ 500ppm
Sulfate	SO ₄ ²⁻	No interference for SO ₄ ²⁻ ≤ 500ppm
Sulfite	SO ₃ ²⁻	Interference for SO ₃ ≥ 20 ppm, giving COD high by 90% at 250 ppm SO ₃ ²⁻
Sulfide	S ²⁻	Interference for S ²⁻ > 0 ppm, giving COD high by >100% at 50 ppm S ²⁻

Cations	Formula	Remark
Aluminum	Al ³⁺	No interference for Al ³⁺ ≤ 500ppm
Calcium	Ca ²⁺	No interference for Ca ²⁺ ≤ 500ppm
Chromate	Cr ³⁺	Interference for Cr ³⁺ > 2 ppm, giving low COD
Ferric Iron	Fe ³⁺	No interference for Fe ³⁺ ≤ 500ppm
Ferrous Iron	Fe ²⁺	Interference for Fe ²⁺ > 100 ppm, giving low COD
Magnesium	Mg ²⁺	No interference for Mg ²⁺ ≤ 500ppm
Potassium	K ⁺	No interference for K ⁺ ≤ 500ppm
Silver	Ag ⁺	Interference for Ag ⁺ > 10 ppm, giving low COD
Sodium	Na ⁺	No interference for Na ⁺ ≤ 500ppm
Zinc	Zn ²⁺	No interference for Zn ²⁺ ≤ 500ppm

For further details, please contact the MANTECH technical team at support@mantech-inc.com.