

BOD 5210 Calculations Explained

Final BOD results are reported as the average of all sample dilutions that meet the passing criteria.

1. When results meet criteria, DO final >1.0mg/L and DO depletion > 2.0mg/L:

- $BOD = (D1-D2) - SV_{seed} / P$

D1 = DO of sample on day 1 (mg/L)

D2 = DO of sample on day 5 (mg/L)

S = seed correction = average DO depletion per mL of seed ( $\Delta DO/mL$ )

Vseed = volume of seed in bottle (mL)

P = volumetric fraction of sample used =  $(V_{sample}/V_{bottle}) * PD$

Vsample = volume of sample (mL)

Vbottle = volume of BOD bottle = 300mL (mL)

PD = predilution factor

Example:

Seed correction: 0.32

Sample	D1	D2	Depletion	Sample Vol	Seed Vol	PD	BOD	BOD Average
GGA_BOD	5.71	4.4	1.31	2	1	1	148.5	
GGA_BOD	7.29	5.1	2.19	4	1	1	140.25	
GGA_BOD	6.58	3.5	3.08	6	1	1	138	
GGA_BOD	6.83	2.7	4.13	8	1	1	142.875	140.37

2. When results fail final DO criteria, i.e. DO2 <1.0mg/L:

- Select the bottle with the highest DO2 value and report as:
- $BOD > (D1-D2) - SV_{seed} / P$

Example:

Seed correction: 0.32

Sample	D1	D2	Depletion	Sample Vol	Seed Vol	PD	BOD	BOD Average
GGA_BOD	5.23	0.3	4.93	2	1	1	691.5	
GGA_BOD	5.11	0.4	4.71	4	1	1	329.25	
GGA_BOD	5.13	0.6	4.53	6	1	1	210.5	
GGA_BOD	5.14	0.9	4.24	8	1	1	147	>147

3. When results fail DO depletion criteria, i.e.  $D1-D2 < 2.0\text{mg/L}$ :

- Select the bottle with the least diluted sample.
- Typically, the bottle with the highest sample volume but watch out for pre-dilution values!
- Select the bottle with the least diluted sample volume, considering pre-dilution.
- $BOD < (D1-D2) - SV_{\text{seed}} / P$

Example:

Seed correction = 0.32

Sample	D1	D2	Depletion	Sample Vol	Seed Vol	PD	BOD	BOD Average
GGA_BOD	5.71	4.8	0.91	2	1	1	88.5	
GGA_BOD	7.29	6.43	0.86	4	1	1	40.5	
GGA_BOD	6.58	6.04	0.54	6	1	1	11	<11
GGA_BOD	6.83	5.84	0.99	8	1	2	50.25	

Least diluted sample = largest number when sample volume/PD

GGA\_BOD 1 =  $2/1 = 2$

GGA\_BOD 2 =  $4/1 = 4$

GGA\_BOD 3 =  $6/1 = 6$

GGA\_BOD 4 =  $8/2 = 4$

Therefore, use GGA\_BOD 3 to calculate BOD:

$BOD < (6.58-6.04)-(0.32*1) / (6/300) *1$

$BOD < 11\text{mg/L}$

**All BOD average results that use failing dilutions, i.e. conditions 2. and 3. must be flagged in the report.**