

# TECHNICAL BULLETIN


## NUMBER 2021 - 012

**Date:** November 2, 2022  
**Author:** James Howley  
**Subject:** Resetting MiniHub to factory defaults

If a MiniHub will not establish communication using the default IP address of 192.168.1.10, a reset to factory defaults may be necessary to resolve this issue.

To do the reset, follow the instructions below:

1. Turn off the power to the MiniHub.
2. Undo the 4 screws securing the cover and remove the cover.
3. Locate switch SW2 on the board near the network jack and toggle the switch to the other position. It may be stiff, so a small screwdriver may help with toggling it.



4. Turn the power back on to the MiniHub. You should see alternating green and red lights, LED 2 and LED3.
5. Wait 10 seconds and then turn off the power again.
6. Toggle the switch back to its original position and turn the power back on.
7. The green light, LED3, should blink on and off once per second.

8. Try pinging 192.168.1.10 or establish communication using MANTECH software.
9. If pinging worked or communication was established, the MiniHub is functional again.
10. Resecure the top cover.

If this does not resolve the issue, please contact the MANTECH technical team at [support@mantech-inc.com](mailto:support@mantech-inc.com).