

The SmartCal™ system is the smart way to automate. Calibrations are time consuming and occupy space that could be used for samples. With SmartCal, calibrations are pre-scheduled, can be run unattended and are completely automated.

The SmartCal Module:

- **Primes all pumps**
- **Rinses the entire Man-Tech System**
- **Calibrates all probes and meters**

At the same time everyday, completely unattended and before chemists arrive in the laboratory!

The SmartCal system will automatically run calibrations at pre-scheduled times. For example, at 6a.m., the Man-Tech system has primed all the pumps, rinsed and the SmartCal system calibrates all electrodes and meters. Upon arrival at the laboratory, you can *immediately* begin running real samples. Calibration and verification results are saved and available for viewing, printing and chemist visual verification at any time.

Key Features and Benefits:

• Pre-scheduled, automated calibrations for:

- pH
- conductivity
- turbidity
- color
- phosphate
- ISE's (fluoride, ammonia, etc.)
- and more.....

• Uses the calibrants from stock bottles

- Ensures consistent calibrations day in and day out
- Calibrants do not occupy autosampler rack space
- Calibrants run operator-free with TitraSip™ SA Systems
- Increases rack capacity by up to 30%

• Improves sample analysis turnaround times

- Begin running real samples immediately in the morning
- No need to wait for the calibrations, they have been completed!

• Use SmartCal to add QC standards during the run

- Increase available rack space for real samples!
- More efficiency
- Up to 9 different calibrants / standards per SmartCal Module
- Up to 4 SmartCal Modules may be connected to a single PC-Titrate™, TitraSip SA, AssayPlus™ or RapidDuo™ System.

• Analyze more real samples

- Increases profit
- Improves cost
- Operate more efficiently

