



## **TECHNICAL BULLETIN**

NUMBER 2018 - 005

Date:February 2018Author:Maggie GriersonSubject:PeCOD® Ranges

The following chart summarizes the COD ranges and sample to electrolyte ratios corresponding to each PeCOD<sup>®</sup> colour range.

PeCOD <sup>®</sup> Range	COD Range (ppm)	Sample: Electrolyte Ratio	Recommended Sample Volume (mL)	Recommend Electrolyte Volume (mL)
Advanced Blue	< 25	3:1	5	15
Green	< 150	1:1	10	10
Yellow	< 1,500	1:9	2	18
Red	< 15,000	1:49	0.5	24.5

## OPTIMIZE YOUR RESULTS. PROTECT OUR ENVIRONMENT.