

TitraMini pH Electrode

NEW Small Diameter pH electrode for MANTECH's MT Series Environmental Titration Systems

Part Number Information

PCE-80-PH1205

TitraMini pH Electrode Glass Double Junction Refillable, 6mm x 100mm, S7 Connector, Full pH Range, No Cable, No Fill Solution. This new small diameter pH electrode is ideal for low sample volumes and skinny tubes. Takes less space in vials than standard 12mm pH electrodes. Well suited for low ionic strength as well as high pH applications. The sample pH is normally registered in less than 10 seconds, and the pH value is reproducible over the entire pH range from 0 to 14 pH.

Electrode Cable

PCE-86-EX1001

Electrode Cable, 6', Detachable



Pictured with TitraPro on left, TitraMini on right