

## TitraPRO4 pH Electrode

*The Superior pH electrode for MANTECH's MT Series Environmental Titration Systems*

### Part Number Information

PCE-80-PH1200D

TitraPRO4 pH Electrode, Comb., DJ, Refillable (KCl), 12mm Dia. X 120mm, S7 Connector, Full pH Range, No Cable. This new Green TitraPRO glass formula eliminates wasted time trying to select the best pH electrode for the job. Green formulation is perfect in low ionic strength solutions as well as high pH applications and is specially formulated to have lightning-fast response time with superior accuracy and repeatability. The sample pH is normally registered in less than 10 seconds, and the pH value is reproducible over the entire pH range from 0 to 14 pH. For the best pH electrode, just think MANTECH green.

### Electrode Cable

PCE-86-EX1001

Electrode Cable, 6', Detachable





Pictured with TitraPro on left, TitraMini on right