

TECHNICAL BULLETIN

NUMBER 2020 - 002

Date: February 26, 2021
Author: Kevin Tularam
Subject: Required Hardware Upgrades for MANTECH Pro Software

MANTECH Pro will only operate on computers with Windows 10 operating system (Home, Professional, or Enterprise).

The following are upgrades that will be required to transition from PC-Titrate to MANTECH Pro:

<i>Type</i>	<i>Required Upgrade</i>
<i>4310 Conductivity Meter</i>	4510 Conductivity Meter + Serial Port
<i>4510 Conductivity Meter</i>	Serial Port Adapter (MANTECH part numbers: 1-Port PC-1250, or 2-Port: PC-1206, depending on system configuration)
<i>Burivar I/II</i>	MANTECH Buret + Ethernet Port
<i>Turbidity AssayPlus Meter</i>	T10 Turbidity Meter
<i>9-Pin Serial Ports</i>	1- or 2-Port USB to Serial adapter (depending on system configuration), must be TrippLite brand, 2-Port adapter must have COM retention (MANTECH part numbers: 1-Port PC-1250, 2-Port: PC-1206)
<i>Ethernet Port</i>	Minimum of 1 – reserved for AutoMax300 Series (or additional USB port for Ethernet adapter). (MANTECH part numbers: USB to Ethernet: PC-1000-1054, D-LINK 5-port (Multiple Burets): PC-1000-1072).
<i>Peripherals</i>	External 4-port SuperSpeed USB 3.0 Hub (MANTECH part number: PC-1000-1128)

Note: Both conductivity meters and Burets do not directly connect to the interface but rather through USB connections. Therefore, a multi-port USB hub may be required, depending on system configuration.

For further details, please contact the MANTECH technical team at support@mantech-inc.com.