

MT3 ANALYZER

OPTIMIZE YOUR RESULTS. PROTECT OUR ENVIRONMENT.



THE MT3 IS THE SIMPLEST SYSTEM IN THE MT SERIES OF ENVIRONMENTAL ANALYZERS. THE MT3 IS DESIGNED FOR NON-TITRATION APPLICATIONS AND NO AUTOMATED SAMPLER, ALTHOUGH IT CAN BE UPGRADED WITH THE ADDITION OF BURETS (TITRATIONS) AND AUTOSAMPLER. THE MT3 IS AN ION METER THAT DOES MORE - INCLUDING PH, CONDUCTIVITY, ISE, TURBIDITY, DO, COD, BOD AND "OPTIMIZED" TOC.

THIS SYSTEM INCLUDES MANTECH PRO™ SOFTWARE AND MANTECH MINIHUB™.

MANTECH Pro™ SOFTWARE

THE NEW SMART SOFTWARE PLATFORM FOR MT ANALYZERS

- INCLUDES BUILT-IN CAPABILITY FOR PH, EC, ISE, TITRATION; TURBIDITY, COLOR, DISSOLVED OXYGEN AND PECOD
 - ALL ENVIRONMENTAL AND SOIL APPLICATIONS EXCEPT BOD5
- DESIGNED FOR WIN10 OS, INCLUDING FULL SECURITY FEATURES FOR NETWORK CONNECTIVITY

MINIHUB™

SMART BI-DIRECTIONAL SOFTWARE-CONTROLLED INTERFACE

- CONNECTS DIRECT TO THE COMPUTER VIA ETHERNET
- SPECIFICATIONS:
 - 4X GALVANICALLY SEPARATED ELECTRODES PORTS WITH BNC CONNECTORS
 - 3X 24V STIRRERS, SOFTWARE CONTROLLED AND VARIABLE SPEED
 - 6X DIGITAL OUTPUTS AND 6X DIGITAL INPUTS
 - 3X LIQUID LEVEL SENSORS FOR BOD