

NEW MT5 AUTOTITRATOR

MAKING WATER QUALITY ANALYSIS SIMPLER

MANTECH PRO

MINIHUB INTERFACE

BURET

STIRRER
PH ELECTRODE (SOLD SEPARATELY)
CONDUCTIVITY PROBE
(SOLD SEPARATELY)

THE NEW MT5 WAS DESIGNED WITH OPERATORS FOR DRINKING WATER AND WASTEWATER TREATMENT PLANTS TO MAKE WATER QUALITY ANALYSIS SIMPLER. THE MT5 IS AN AUTOTITRATOR THAT DOES MORE – WITH THE ABILITY TO ADD MULTIPLE BURETS, CONDUCTIVITY, DO, TURBIDITY, ISE, COLOR, COD, BOD, “OPTIMIZED” TOC AND CAN BE UPGRADED TO AN AUTOSAMPLER.

THIS SYSTEM INCLUDES MANTECH PRO™ SOFTWARE, MANTECH MINIHUB™, BURET, STAND AND STIRRER.

MANTECH Pro™ SOFTWARE

THE NEW SMART SOFTWARE PLATFORM FOR MT ANALYZERS

- INCLUDES BUILT-IN CAPABILITY FOR PH, EC, ISE, TITRATION, TURBIDITY, COLOR, DISSOLVED OXYGEN AND PECOD
- ALL ENVIRONMENTAL AND SOIL APPLICATIONS EXCEPT BOD5
- DESIGNED FOR WIN10 OS, INCLUDING FULL SECURITY FEATURES FOR NETWORK CONNECTIVITY

MINIHUB™

SMART BI-DIRECTIONAL SOFTWARE-CONTROLLED INTERFACE

- CONNECTS DIRECT TO THE COMPUTER VIA ETHERNET
- SPECIFICATIONS:
 - 4X GALVANICALLY SEPARATED ELECTRODES PORTS WITH BNC CONNECTORS
 - 3X 24V STIRRERS, SOFTWARE CONTROLLED AND VARIABLE SPEED
 - 6X DIGITAL OUTPUTS AND 6X DIGITAL INPUTS
 - 3X LIQUID LEVEL SENSORS FOR BOD

SPECIAL
PRICE

ACT FAST!

TAKE ADVANTAGE OF OUR SPECIAL RELEASE PACKAGE PRICE INCLUDING A PH ELECTRODE AND FREE LAPTOP.
EXPIRES DECEMBER 31, 2021

CONTACT US AT:
MANTECH-INC.COM
+1 (519) 763-4245
INFO@MANTECH-INC.COM

 @MANTECH-INC

 LINKEDIN.COM/COMPANY/MANTECH-INC

5473 HIGHWAY 6 NORTH
GUELPH, ON CANADA
N1H 6J2

 **MANTECH**
OPTIMIZE YOUR RESULTS. PROTECT OUR ENVIRONMENT.