

## INTELLIPUMP™

*MANTECH Peristaltic Pump*

**Part Number:** PC-1000-1137

**Product Description:** The MANTECH variable speed peristaltic pump is an extremely versatile pump, driven by a long-lasting, brushless DC motor rated for continuous duty. Featuring IoT connectivity, the pump speaks with MANTECH software to establish connection and control the flow rate and direction.



### Features and Benefits:

- Flow rates of 12.5 to 860 mL/min (depending on tubing diameter)
- 24 VDC
- USB connection
- Accepts continuous tubing for clean fluid paths
- Heavy duty steel rotor
- Friction-reducing Nylatron rollers for longer tube life
- Up to 25 psi working pressure