

MANTECH has relaunched PeCOD®'s *Just Add Water* reagents! Containing all solutes typically present in our previous solutions (now discontinued), customers are only required to add deionized (DI) or distilled (found at your local grocery store) water – aiming to reduce costs and your carbon footprint. To ensure access to all required supplies, MANTECH is offering a special supply kit for existing PeCOD® customers. All related part numbers (e.g., supply kit and solution kits) are listed below.

Please contact orders@mantech-inc.com to place your order today!

EXISTING CUSTOMER SUPPLY KIT

PQK-202-0093		
Part Number	Qty	Description
PC-5000-1941	1	250mL plastic squeeze bottle
PC-5000-1940	1	½ gallon jug
PC-5000-1938	1	500mL graduated cylinder
PC-5000-1948	1	Funnel
PC-5000-1370	4	1L plastic sample bottles with caps
PC-6000-064	2	¼" Grommet (for 1L bottle caps) for waste and QC check standard bottles
PC-5000-1947	1	Straw for bottle top dispenser
		Set of labels for empty bottles (electrolyte, calibrant, blank, QC check standards and waste)

SOLUTION KITS

ADVANCED BLUE RANGE (<20MG/L)

PQK-202-0100			
Part Number	Qty	Volume	Description
PQA-201-0050	1	1L	Advanced Blue Range Electrolyte
PQA-201-0070	1	250mL	Advanced Blue Range Blank
PQA-201-0060	1	250mL	Advanced Blue Range Calibrant
PQA-201-0080	1	250mL	Advanced Blue Range QC Check

SOLUTION KITS (CONT.)

GREEN RANGE (<150MG/L)

PQK-202-0101			
Part Number	Qty	Volume	Description
PQA-201-0051	1	1L	Green Electrolyte
PQA-201-0071	1	250mL	Green Blank
PQA-201-0061	1	250mL	Green Calibrant
PQA-201-0081	1	250mL	Green QC Check

YELLOW RANGE (<1,500MG/L)

PQK-202-0102			
Part Number	Qty	Volume	Description
PQA-201-0052	1	1L	Yellow Electrolyte
PQA-201-0072	1	250mL	Yellow Blank
PQA-201-0062	1	250mL	Yellow Calibrant
PQA-201-0082	1	250mL	Yellow QC Check

RED RANGE (<15,000MG/L)

PQK-202-0103			
Part Number	Qty	Volume	Description
PQA-201-0052	1	1L	Red Electrolyte
PQA-201-0073	1	250mL	Red Blank
PQA-201-0063	1	250mL	Red Calibrant
PQA-201-0083	1	250mL	Red QC Check

Please contact orders@mantech-inc.com to place your order today!