


## INTRODUCING PECOD'S 'JUST ADD WATER' SOLUTIONS

Containing all solutes typically present in our previous solutions, the *Just Add Water* solutions are concentrated into smaller volumes requiring the addition of water (deionized or distilled). By adding water after shipping, we have cut the volume and weight of the bottles by up to 75% thereby reducing plastic waste, shipping costs, and carbon emissions - ***ultimately reducing your cost per sample with a green impact!***



REDUCED PLASTIC WASTE



CUT TRANSPORT COSTS



IMPROVED PRICING



EASY SAMPLE PREP

### DO YOU HAVE QUESTIONS OR WANT MORE INFORMATION?

Visit our website for:

- *Easy-to-follow* solution preparation instructions and tutorial videos
- *Step-by-step* sample preparation tutorial videos
- Frequently asked questions
- and more!

**VISIT [MANTECH-INC.COM](http://MANTECH-INC.COM)**

# SOLUTION KITS

For each range, solution kits will include 250mL blank slurry, 250mL calibrant slurry and 250mL QC check slurry, 500mL (or 1L for Adv. Blue range) electrolyte slurry. The blank, calibrant, and QC check slurry will make up 1L worth of solution following the addition of the deionized or distilled water. The electrolyte slurry will make up 2L worth of solution following the addition of deionized or distilled water.

*All solutions can be purchased separately or in a kit.*

## ADVANCED BLUE RANGE (COD < 20mg/L)

PQK-202-0100				
Part Number	Qty	Description	Shipped Volume	Total Volume (w/ water)
PQA-201-0050	1	Adv. Blue Electrolyte	1L	2L
PQA-201-0070	1	Adv. Blue Blank	250mL	1L
PQA-201-0060	1	Adv. Blue Calibrant	250mL	1L
PQA-201-0080	1	Adv. Blue QC Check	250mL	1L

## GREEN RANGE (COD < 150mg/L)

PQK-202-0101				
Part Number	Qty	Description	Shipped Volume	Total Volume (w/ water)
PQA-201-0051	1	Green Electrolyte	500mL	2L
PQA-201-0071	1	Green Blank	250mL	1L
PQA-201-0061	1	Green Calibrant	250mL	1L
PQA-201-0081	1	Green QC Check	250mL	1L

## YELLOW RANGE (COD < 1,500mg/L)

PQK-202-0102				
Part Number	Qty	Description	Shipped Volume	Total Volume (w/ water)
PQA-201-0052	1	Yellow Electrolyte	500mL	2L
PQA-201-0072	1	Yellow Blank	250mL	1L
PQA-201-0062	1	Yellow Calibrant	250mL	1L
PQA-201-0082	1	Yellow QC Check	250mL	1L

## RED RANGE (COD < 15,000mg/L)

PQK-202-0103				
Part Number	Qty	Description	Shipped Volume	Total Volume (w/ water)
PQA-201-0052	1	Red Electrolyte	500mL	2L
PQA-201-0073	1	Red Blank	250mL	1L
PQA-201-0063	1	Red Calibrant	250mL	1L
PQA-201-0083	1	Red QC Check	250mL	1L

# SUPPLY KIT

## SPECIAL PACKAGE FOR EXISTING CUSTOMERS

To ensure access to all required supplies, MANTECH is offering a special supply kit for existing PeCOD® customers. Includes 2L bottle (electrolyte), straw for bottle top dispenser, funnel, graduated cylinder, squeeze bottle, 4x 1L bottles (calibrant, blank, QC check, waste). All related part numbers are listed below.

*Supply kits do not include reagents. Add selected range(s) of reagents separately.*

### Supply Kit for Existing PeCOD® Customers PQK-202-0093

#	Part Number	Qty	Description
1	PC-5000-1941	1	Plastic squeeze bottle
2	PC-5000-1370	4	1L plastic sample bottles
3	PC-5000-1940	1	2L jug (for electrolyte)
4	PC-5000-1948	1	Funnel
5	PC-5000-1947	1	Straw for bottle top dispenser
6	PC-5000-064	2	1/4" grommet (for 1L bottle caps) for waste and QC check standard bottles
7	PC-5000-1938	1	500mL graduated cylinder
			Set of labels for empty bottles (electrolyte, blank, QC check standards, and waste)



Have Questions or Want More Information?

[MORE INFORMATION](#)

Need to Place an Order?

[ORDER HERE](#)



[WWW.MANTECH-INC.COM](http://WWW.MANTECH-INC.COM)

[INFO@MANTECH-INC.COM](mailto:INFO@MANTECH-INC.COM)

+1 (519) 763-4245

+1 (866) 763-2122

[@mantech\\_inc](https://twitter.com/mantech_inc)

[linkedin.com/company/mantech-inc](https://www.linkedin.com/company/mantech-inc)

© 2022 MANTECH. All rights reserved.

PeCOD® is a registered trademark of  
MANTECH.