


# NEW AM430 Series Autosamplers

For use with MANTECH's Automated MT Series & BOD Pro Systems



- ✓ Utilizes Robust Robotics
- ✓ Maximizes Limited Bench Space
- ✓ Automates Low Throughput Analyses


**Fits comfortably on lab benches of 24 inch depth!**

## Standard AM430 Rack Capacities

*Custom racks available upon request.*

|  |  |  |  |
|--|--|--|--|
| <b>15mL Tubes</b><br>99 to 495 positions | <b>50mL Tubes</b><br>68 to 340 positions | <b>125mL Cups</b><br>20 to 100 positions | <b>300mL Beakers</b><br>12 to 60 positions |
|--|--|--|--|




## AM430 Series Autosamplers Dimensions



*\*Ideal for labs with space restrictions, the AM430s save 5 inches in depth compared to the AM400 autosamplers.*



Consult our team about the autosampler that suits your requirements best!

 [www.mantech-inc.com](http://www.mantech-inc.com)  
 519-763-4245  
 +1 866-763-4245 (toll free)  
 [info@mantech-inc.com](mailto:info@mantech-inc.com)