



MANTECH

OPTIMIZE YOUR RESULTS. PROTECT OUR ENVIRONMENT.



5-DAY BIOCHEMICAL OXYGEN DEMAND (BOD5) SOFTWARE PACKAGE



COMPLIANCE MEETS FLEXIBILITY

BOD Pro™ is lab automation software coded and designed for MANTECH analyzers. Offered with automation or as a stand-alone software, BOD Pro™ is built with features to simplify and streamline 5-day biochemical oxygen demand (BOD5) analyses while delivering valid and compliant results.

SOFTWARE ONLY PACKAGE
\$2,410 USD



BOD PRO™ CONFORMS TO THE FOLLOWING REGULATORY NORMS:

- ASTM STANDARD METHODS 5210B, 23RD AND 24TH EDITIONS
- EPA 405.1
- ISO EN 5815-1
- ISO EN 1899-1, -2

Probe Options

Membrane or optical

Single, dual or triple probes

YSI MultiLab 4010-W integration or manual entry required

SOFTWARE FEATURES



Customizable Shortcut Buttons

- Simplifies and reduces prep on routine sample runs, maintenance setup and more



IoT Device Connectivity

- Bi-directional communication with hardware and software



Import & Export LIMS Capabilities

- Import sample run information
- Auto-exporting feature to LIMS or network location



Enable OS & Security Updates

- Designed for Windows 10 and 11 Pro OS

UPGRADE TO BOD PRO™ TODAY!

Plan on upgrading to Windows 10 or 11 Pro? Due to Microsoft's major security updates, PC-BOD is not guaranteed to be operational with the latest OS updates. MANTECH's **NEWEST** BOD Pro™ software is designed for Windows 10 & 11 Pro OS for enhanced cyber security and dependability.



OPTIMIZED SECURITY

- Flexible user account system with easy configuration and user defined permissions



ENHANCED USER OPERATION

- Intuitively designed with colour coding, visual rack display and live script monitoring features



HARDWARE INTEGRATION

- PC-BOD hardware is not compatible with BOD Pro™ software and requires MANTECH's [LATEST hardware](#)

CONFIGURATIONS TO OPTIMIZE YOUR THROUGHPUT

SIMPLIFY MANUAL DO READINGS




- ✓ Automated calculations
- ✓ Incubation period tracking
- ✓ Use your existing DO meter

AUTOMATE DO READINGS AND SAMPLE PREPARATION




- ✓ Up to 108 position autosamplers
- ✓ Optional dilution water, seed and/or inhibitor addition


BOD PRO™ SOFTWARE




Home Screen With Customizable Shortcut Buttons



Up to Two Single, Dual or Triple Probe DO Meters



Build Run Templates With Intuitive Rack Display



Add Up to 10 User-Defined Columns in Report Templates



5473 HIGHWAY 6 NORTH
GUELPH, ONTARIO, CANADA
N1H 6J2

WWW.MANTECH-INC.COM
INFO@MANTECH-INC.COM
+1 (519) 763-4245
+1 (866) 763-2122

 @mantech_inc
 linkedin.com/company/mantech-inc-

© 2025 MANTECH. All rights reserved.
PeCOD® is a registered trademark of
MANTECH.